

INspire UT



Our **INspire™ Ultra-Translucent** Multi-Layer Zirconia is ideal for both anterior and posterior units, including up to 3-unit anterior or posterior bridges, and in some cases veneers. It is fabricated to provide smooth shade transitions, which mimic the appearance of natural tooth enamel, dentin, and cervical shades. Since the color progression is natural and undisrupted, there are no abrupt shade edges. Tooth-like preshading provides an efficient and economical process that eliminates the need for any staining procedures. It has a high bending strength of up to 850 MPa, which is comparable to alternative all-ceramic restorations. This zirconia can match all 16 VITA® Classical Shades. For optimal esthetics, please provide the stump shade with each INspire™ UT case.

Indications

INspire UT can be used for almost any C & B or screw-retained implant case when anterior or posterior crowns, veneers and 3-unit bridges are prescribed. Also it can be used as an esthetic alternative to a PFM with metal occlusion due to limited occlusal or lingual space.

Contraindications

When esthetic expectations are high and it is important that the restorations match surrounding natural dentition or other existing restorations. If bonding is necessary to retain the restoration, bond strength is weaker and less predictable than other ceramics.

Prep Guides

Shoulder preparation not needed. A mild chamfer or a feather-edge margin is good. 1mm buccal, lingual, and occlusal reduction is ideal, but can go to .5mm in some areas when reduction is limited. Minimum occlusal reduction of 0.5 mm; 1 mm is ideal. Adjustments and polishing: Adjust full-contour zirconia crowns and bridges using water and air spray to keep the restoration cool and to avoid microfractures with a fine grit diamond. If using air only, use the lightest touch possible when making adjustments. A football-shaped bur is the most effective for occlusal and lingual surfaces (on anterior teeth); a tapered bur is the ideal choice for buccal and lingual surfaces. Polish full-contour zirconia restorations with the porcelain polishing system of your choice.

Cementation

Full-contour zirconia may be cemented using a resin reinforced glass ionomer such as RelyX Luting cement. When a greater bond is needed do to the lack of a retentive preparation, use resin cement like RelyX Unicam or RelyX Ultimate. Before cementing all full-contour zirconia crowns, the interior surface of the crown needs to be cleaned with Ivoclean (Ivoclar Vivadent - Amherst, NY). This is critical in assuring maximum bond strength.