

For Each Indication — The Ideal Cement



RelyX™ Ultimate
Adhesive Resin Cement



RelyX™ Unicem
RelyX™ Unicem 2
Self-Adhesive Resin Cement



RelyX™ Luting Plus
Resin Modified Glass
Ionomer Cement



RelyX™ Veneer
Veneer Cement



	RelyX™ Ultimate	RelyX™ Unicem	RelyX™ Luting Plus	RelyX™ Veneer
Metal / Metal based				
Inlays / Onlays	+	++	+	-
Crowns / Bridges	+	++	++	-
Endodontic Posts	+	++	++	-
Maryland Bridges	++	+	-	-
On Implant Abutments	+	++	++	-
Glass Ceramics (incl. Li. Disilicate) (e.g. e.max, VITA Mark II, IPS Empress 2)				
Inlays / Onlays / Table Tops	++	+	-	-
Crowns / Bridges	++	++	-	-
Maryland Bridges	++	+	-	-
Veneers	+	-	-	++
On Implant Abutments	+	+	-	-
Oxide Ceramics (e.g. Lava™ Zirconia, Brux Zir, Procera)				
Inlays / Onlays	++	++	+	-
Crowns / Bridges	+	++	+	-
Endodontic Posts	+	++	+	-
Maryland Bridges	++	+	-	-
On Implant Abutments	+	++	++	-
Resin Nano Ceramics (e.g. Lava™ Ultimate CAD/CAM Restorative)				
Inlays / Onlays	++	-	-	-
Crowns	++	-	-	-
Veneers	+	-	-	++
On Implant Abutments	+	-	-	-
Resin Composites				
Inlays / Onlays	++	+	-	-
Crowns	++	++	-	-
Endodontic Posts	+	++	-	-
Veneers	+	-	-	++
On Implant Abutments	++	++	-	-

- ++ Highly recommended*
- + Recommended
- Not indicated


Delivery Method


Automix Syringe


Clicker™ Dispenser


Capsule


Powder/Liquid


Syringe

* Either better performance for this indication or easier handling with equal performance.

Pretreatment of Restorations



Glass Ceramic: E-max, IPS Empress 2, Vita Mark II,
Veneers

Pretreat with Hydrophloric Acid

Add Silane to surface with Scotchbond Universal,
Monobond, etc

(Optional: Ivoclean to clean Intaglio Surface of the crown)

Oxide Ceramic: Lava Plus, Brux Zir, Procera

Pretreat with Air Abrasion using Aluminum Oxide

Add Silane to surface with Scotchbond Universal,
Monobond, etc

(Optional: Ivoclean or Air Abrasion to clean Intaglio Surface
of the crown)