



Product

Comparison Guide



*Most Popular

| | INspire HS High Strength Zirconia | INspire UT Ultra Translucent Zirconia | IPS e.max |
|------------------------|---|---|------------------------------|
| IDEAL SITUATION | Posterior | Anterior or posterior | Single-unit crown or veneer |
| | Cases requiring additional durability, including those for the molar region | Cases requiring additional translucency, layered look, and esthetics over Inspire HS, but with almost equivalent strength | Patient in optimal occlusion |
| INDICATIONS | Crowns | Crowns | Crowns |
| | Bridges: 6+ unit to full-arch posterior bridges | Bridges: 6+unit full-arch anterior and posterior bridges | |
| | Inlays / Onlays | Inlays / Onlays | Inlays / Onlays |
| | | Veneers | Veneers |
| | Implants: Singles, bridges with any combination of abutments and pontics, inlay bridges, Maryland type bridges and screw-retained implants | Implants: Singles, bridges with any combination of abutments and pontics, inlay bridges, Maryland type bridges and screw-retained implants | |
| TRANSLUCENCY | ★★★ | ★★★★★ | ★★★★★ |
| STRENGTH | 1300 MPa | 1200 MPa | 400 MPa |
| DURABILITY | ★★★★★ | ★★★★★ | ★★★ |
| MATERIAL | Zirconium dioxide | Zirconium dioxide | Lithium disilicate |



R-Dent
DENTAL LABORATORY



901.372.8020



www.rdentlab.com



7490 Bartlett Corporate Cove W, Bartlett, TN 38133