

BECAUSE FIRST IMPRESSIONS COUNT

Bite Registration

Issue

- Difficulty in establishing vertical and horizontal relationships of the dental arches to the mounted stone casts.
- Excessive occlusal adjustment.

Cause

- Bite registration not provided.
- Poor interocclusal record.

Solution

- Provide a rigid and dimensionally stable bite registration with each case.
- Properly trim a good bite registration material.
- Use dimensionally stable bite registration material.
- Ensure that the patient does not move during the procedure.

For additional assistance contact the 3M ESPE Technical Hotline: 1-800-634-2249 In Canada, call 1-800-265-1800 x6229