



Zirlux Zirconia *Cementation and Bonding Reference*

Zirlux® restorations may be cemented or bonded using conventional cements and luting materials. Resin modified glass ionomers and resin cements are recommended over zinc oxide phosphate cements to optimize aesthetic results. Inlays, onlays, and veneers must be adhesively cemented (bonded).

Ceramir® Crown & Bridge

Doxa

RelyX™ Luting Cement

3M ESPE

GC Fuji Plus®

GC

Panavia™ F2.0

kuraray

Multilink® Automix

**ivoclar
vivadent**
passion vision innovation



Refer to the respective manufacturer's instructions for use and proper application of any cementation/bonding materials.



Discover why Zirlux is *right for you.*

www.zirlux.com

#iamtypez



7490 Bartlett Corporate Cove
West Bartlett, TN 38133

877-RDENT4U

www.rdentlab.com